



ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

Corporate Headquarters
AMRESKO, LLC
6681 Cochran Road
Solon, Ohio USA 44139-4300
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

CERTIFICATE OF QUALITY

Product: BIURET PROTEIN ASSAY REAGENT
Grade: PROTEOMICS GRADE
Code: M262

Appearance: Blue clear liquid.

Mol. Formula: Not Applicable

Mol. Weight: Not Applicable

Description:

<u>Assay:</u>	<u>Min</u>	<u>Max</u>	<u>Units</u>	<u>Disposition</u>
Abs. @ Lambda Max	0.31	0.37		PASS
Lambda Max	672	676	nm	PASS

Storage: Store at room temperature. Protect from moisture.

M262PTM

This signature indicates that the above material has met all quality specifications and has been reviewed by a quality representative.

Signature:

Title:

QC MANAGER

Date:

6/4/2012