
Cat.-No.	Size
AB-105S	1 ml
AB-105L	5 ml

Description:

Contains 50 mg/ml LEXSY Neo antibiotic.
Sterile. Based on antibiotic Geneticin™.

Storage conditions:

Store at -20°C.
Stable for 18 months.

Selection:

For selection of recombinant *Leishmania* strains
LEXSY Neo is added to growth medium to a
final concentration of 50 µg/ml.

Find more information in the User's Guides for
LEXSY kits (download pdf from the LEXSY section
of our website).

FOR RESEARCH USE ONLY.

**NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC
OR THERAPEUTIC USES.**