


FS-IU638A



Sterile Water, Pyrogen Free, Endotoxin Free

Product Information

Name: Sterile Water, Endotoxin Free (apyrogen)

Description: Autoclaved sterile ultrapure water, tested devoid endotoxins.

Cat. Number: IU638A, 1L
IU638B, 4 L

CAS: 7732-18-5

Molecular Formula: H₂O

Molecular Weight: 18.02

Recommended storage: Room Temperature

Sterile Water is reverse osmosis purified to 18Mohms (Resistivity), Bacterial filtered (0.2 micron filtered) and autoclaved sterile. Endotoxin removal is facilitated by electrostatic filtration and the product is LAL tested for pyrogens. The LAL tested water conforms to the standard having less than <0.25EU/ml to ensure the water is of pre-requisite quality. This product is therefore pyrogen free. CFU>0 WFi compatible.

This water is ideal for cell culture applications, and suits most other applications, i.e. molecular biology,...

For any information, please contact / Pour toute information complémentaire, s'adresser à :
Uptima -Interchim, Division BioSciences interbiotech@interchim.com; Hotline : +33 (0)4 70 03 73 06

Rev.U07E-S07