

Product Information

Mini Cell Scrapers


Catalog Number: 22003

Storage: Room temperature

Quantity: 200 scrapers total, 10 scrapers per sterile packet

Product Description

Sterile, polyethylene, disposable mini cell scrapers can be used to effectively harvest cells or cell lysates from 24, 48 and 96 well plates. They measure 3/16" (0.5 cm) wide and 2 3/8" (6 cm) long. Mini cell scrapers may be used with the AccuEasy™ Flow Cytometry Kit (catalog number 30069) for lifting highly adherent cells.



Related Products

Catalog number	Product
30069	AccuEasy™ Flow Cytometry Kit
23006	Flow Cytometry Fixation/Permeabilization Kit
22015	Fixation Buffer
22016	Permeabilization Buffer
22017	Permeabilization and Blocking Buffer
22010	10% Fish Gelatin Blocking Buffer
22011	Fish Gelatin Powder
22014	30% Bovine Serum Albumin Solution
22012	Dry Milk Powder
22002	Tween®-20

Please visit www.biotium.com to view our full selection of products featuring bright and photostable fluorescent CF™ dyes, including secondary antibodies and Mix-n-Stain™ antibody labeling kits. Biotium also offers a variety of apoptosis and cell viability assays for flow cytometry analysis, including mitochondrial membrane potential dyes, fluorescent Annexin V conjugates, and NucView™488 Caspase-3 Substrate for live cells.

Materials from Biotium are sold for research use only, and are not intended for food, drug, household, or cosmetic use.