


FT-JO3540



## Sequencing Gel 19:1

*Sequencing Gel 19:1 ratio*

### Product Description

**Name :** **6% Sequencing Gel 19:1 ratio**  
 6% Manual Sequencing Gel Solution Bis-Acrylamide Ratio 19:1 (Ready to use formulation)  
**Catalog Number :** JO3540, 500ml  
 JO3541, 1L

**Name :** **8% Sequencing Gel 19:1 ratio**  
 8% Manual Sequencing Gel Solution Bis-Acrylamide Ratio 19:1 (Ready to use formulation)  
**Catalog Number :** 1Q6690, 1L  
 1Q6691, 500 ml

**Name :** **4.5% Sequencing Gel 19:1 ratio**  
 4.5% Manual Sequencing Gel Solution Bis-Acrylamide Ratio 19:1 (Ready to use formulation)  
**Catalog Number :** 1R0360, 500 ml  
 1R0361, 1 L

**Storage:** +2-8°C Minimum shelf life 6 months

Used in routine DNA sequencing work, ready to use formulation contains 7M Urea (Denaturant) and 1X TBE buffer. Each lot is DNase/ RNase tested. All stock solutions are made from ultra-pure components and both 0.45 and 0.2micron filtered.

### Ordering information

For any information, please ask : Uptima / Interchim; Hotline : +33(0)4 70 03 73 06

[Order on-line](#) or [Contact](#) your local distributor

**Disclaimer :** Materials from Uptima are sold **for research use only**, and are not intended for food, drug, household, or cosmetic uses. Uptima is not liable for any damage resulting from handling or contact with this product.