

ReadiUse™ Mammalian Cell Lysis Buffer *5X*

Ordering Information:

Product Number: 20012 (10 mL)

Storage Conditions:

Store at 2-8 °C

Expiration date is 6 months from the date of receipt

Typical Assay Protocol

1. Warm up ReadiUse™ mammalian cell lysis buffer *5X* at room temperature before use.
2. Prepare 1X ReadiUse™ mammalian cell lysis buffer by adding **1 ml** of 5X Screen ReadiUse™ mammalian cell lysis buffer into **4 ml** of Mili-Q water, and mix them well.
3. The 1X ReadiUse™ mammalian cell lysis buffer (from Step 2) now is ready to use.
4. Treat cells as desired, wash plate with PBS to remove residual media.
5. The following table is the suggested volume for lysis adherent cells.

Plate	384-well	96-well	48-well	24-well	12-well	6-well	100 cm ²
Volume	20 µL	50 µL	100 µL	150 µL	200 µL	300 µL	800 µL

6. Lyse non-adherent cells by adding about 100 µL of 1X ReadiUse™ mammalian cell lysis buffer to 10⁶ cell pellet.
7. Incubate at room temperature for 5 min.
8. The cell lysate now is ready for further biochemical analysis.

Disclaimer: This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact our technical service representative for more information.

www.interchim.com



211 bis Avenue Kennedy - BP 1140
03103 Montluçon - France
33 (0) 4 70 03 88 55
Fax 33 (0) 4 70 03 82 60
e-mail interchim@interchim.com

Agence Paris - Normandie
33 (0) 1 41 32 34 40
Fax 33 (0) 1 47 91 23 90
e-mail interchim.paris@interchim.com