


FT-LO8990



interchim®

INTERFINE CHEMICALS ANALYTICAL SCIENCES BIOCHROMATOGRAPHY LIFE SCIENCES

+33 4 70 03 73 06 interbiotech@interchim.com

## Thimerosal, 1% w/v solution

### Product Description

**Name :** Thimerosal, 1% w/v solution  
**Catalog Number :** LO8990, 100 ml  
**Structure :** CAS [54-64-8]

**Storage:** Room temperature

Thimerosal solution is widely used as a preservative and deemed to be less hazardous when used as a solution. Suitable for Gel Filtration Chromatography , Affinity Chromatography, Anti-body preservation , Molecular biology, Cell Biology and General Biological use.

1% Thimerosal Solution is generally used and diluted to 0.1% strength with deionised RO purified water, it is also 0.2micron filtered. DNase RNase free.

### Ordering information

Catalog size quantities and prices may be found at <http://www.interchim.com>.  
Please inquire for higher quantities (availability, shipment conditions).

For any information, please ask : Uptima / Interchim; Hotline : +33(0)4 70 03 73 06

[Order on-line](#) or [Contact](#) your local distributor

**Disclaimer :** Materials from Uptima are sold **for research use only**, and are not intended for food, drug, household, or cosmetic uses.  
Uptima is not liable for any damage resulting from handling or contact with this product.

LIFE SCIENCES

**interchim**

211 bis, avenue JF Kennedy  
BP 1140 - 03100 Montluçon  
Fax +33 4 70 03 82 60

■ Hotline +33 4 70 03 73 06 • [interbiotech@interchim.com](mailto:interbiotech@interchim.com) ■