



ISO 9001 CERTIFIED
ISO 13485 CERTIFIED

Corporate Headquarters
AMRESKO, LLC
6681 Cochran Road
Solon, Ohio USA 44139-4300
440/349-1199 FAX: 440/349-1182
www.amresco-inc.com
Email: info@amresco-inc.com

CERTIFICATE OF QUALITY

Product: AGAR, BACTERIOLOGICAL
Grade: BACTERIOLOGICAL GRADE
Code: J637

Appearance: White to light cream powder

Mol. Formula: Not Applicable

Mol. Weight: Not Applicable

Description:

<u>Assay:</u>	<u>Min</u>	<u>Max</u>	<u>Units</u>	<u>Disposition</u>
Ash		6.5	%	PASS
Gel Point (1.5%)	32	38	Degrees C	PASS
Gel Strength (1.5%)	550	950	g/cm2	PASS
Heavy Metals		0.002	%	PASS
Loss on Drying		12.0	%	PASS
Melting Point	80	90	Degrees C	PASS
Pass Thru 60 Mesh	90		%	PASS
pH (1.5%)	6.0	7.5		PASS
Solubility (Boiling Water)	PASS		P/F	PASS
Turbidity		40	Nephelos	PASS

Storage: Store at room temperature. Protect from moisture.

J637BACT

This signature indicates that the above material has met all quality specifications and has been reviewed by a quality representative.

Signature:

Title:

QC MANAGER

Date:

7/10/2012