
Specification Sheet

HEPES Buffer

Product Information

Name:	HEPES Buffer Solution, 1M, pH 7.3, 0.2 μ filtered
Cat. Number:	N1466A, 20 ml
Molecular Formula:	C ₈ H ₁₈ N ₂ O ₄ S
Molecular Weight :	238.30
Recommended storage:	Room Temperature.

Quality control

<i>Test</i>	<i>Specifications</i>
Appearance:	Colorless liquid
pH at 37°C	7.2-7.4
DNase, RNase, Protease activity	None

For any information, please contact / Pour toute information complémentaire, s'adresser à :
Uptima -Interchim, Division BioSciences Interbiotech@interchim.com; Hotline : +33 (0)4 70 03 73 06

Contact your local distributor Uptima, powered by

uptima@interchim.com



213 Avenue J.F. Kennedy - BP 1140
03103 Montluçon Cedex - France
Tél. 04 70 03 88 55 - Fax 04 70 03 82 60