

FT-N14872



M63 Medium Broth

Medium with low osmolarity resulting growth slower cells

Product Description

Name :	M63 Broth
Catalog Number :	N14873, 500 g
Storage:	at +2 to +4°C ^(L) Protect from light and moisture Once opened keep powdered medium closed to avoid hydration.

USES

Suspend 20.012 grams of medium in one liter of distilled water. Mix well and dissolve by heating with frequent agitation. Boil for one minute until complete dissolution. Dispense into appropriate containers and sterilize in autoclave at 121°C for 15 minutes. The medium when cooled is clear.

FORMULA (reconstituted)

Monopotassium phosphate : 3 g/l	Dextrose : 8 g/l	Ammonium Sulfate : 2 g/l
Dipotassium Phosphate : 7 g/l	Magnesium Sulfate : 0.012 g/l	

Related / associated products and documents

- [Cell Culture Media Components](#) ^[PL]
e.g. Agar #291901, Terrific Broth #82111A, NZCYM Broth # N1473b, Marine Broth #-A2WT30, Brain Heart Infusion Broth JI0611, TCBS Agar #CJ2382, ...

- Other cell culture and assay reagents

[FT-N68091](#) Accumax, cells clumps dissociation reagent

[PH-BE007a](#) CosiGel & CosiMatrix 3D Cell Culture

See [Product highlights](#), [BioSciences Innovations catalogue](#) and [e-search tool](#).

Ordering information

For any information, please ask at uptima@interchim.com or at Uptima / Interchim; Hotline : +33(0)4 70 03 73 06

Disclaimer : Materials from Uptima are sold **for research use only**, and are not intended for food, drug, household, or cosmetic uses. Uptima is not liable for any damage resulting from handling or contact with this product. Rev.T08E-T04V