


FT-NJI661



Hanks' Balanced Salt Solution (HBSS)

A sterily filtered saline solution suitable for a variety of cell culture applications

Product Description

Catalog #: NJI661, 500 ml
Name : **Hanks' Balanced Salt Solution (HBSS)**
Lacking phenol red
Use: Dilute to 1X concentration using sterile water.
Storage : 4° C

For Research Use Only

Technical and Scientific Information

HBSS, 10X is a concentrated formula of Hank's Balanced Salt Solution. HBSS can be used in an assortment of cell culture applications, such as diluting cells for counting, preparing reagents, transporting cells or tissues, as well as washing cells before dissociation. HBSS provides a buffering system to maintain the physiological pH range and osmotic balance of the culturing media. This buffer also provides a source of water, essential inorganic ions, and carbohydrates for the cell or tissue culture.

Ordering information

Catalog size quantities and prices may be found at <http://www.interchim.com>.
Please inquire for higher quantities (availability, shipment conditions).

Please contact InterBioTech – Interchim for any other information
Hotline : +33(0)4 70 03 73 06 – Interbiotech@interchim.com

Disclaimer : Materials are sold **for research use only**, and are not intended for food, drug, household, or cosmetic use. Interchim® is not liable for any damage resulting from handling or contact with this product.