


FT-NJK63A



Sodium Azide, 10% solution

Product Description

Name :	Sodium Azide, 10% solution STERILE	
Catalog Number :	NJK63A, 100 ml	NJK631, 1 L
Structure :	CAS: 26628-22-8	
Storage:	Room temperature (Z)	



Sodium Azide solution is widely used as a preservative and deemed to be less hazardous when used as a solution. Suitable for Gel Filtration Chromatography, Affinity Chromatography, Anti-body preservation, Molecular biology, Cell Biology and General Biological use.

10% w/v Sodium Azide Solution is generally used and diluted to 0.01% w/v strength with deionised RO purified water, it is also 0.2micron filtered.

Ordering information

Catalog size quantities and prices may be found at <http://www.interchim.com>.

For any information, please ask : Uptima / Interchim; Hotline : +33(0)4 70 03 73 06

[Order on-line](#) or [Contact](#) your local distributor

Disclaimer : Materials from Uptima are sold **for research use only**, and are not intended for food, drug, household, or cosmetic uses. Uptima is not liable for any damage resulting from handling or contact with this product. T04E

