

FT-RP2327



1,3-Dioleoyl-2-palmitoylglycerol

An triglyceride used in making human milk replacements

Product Description

Name: 1,3-Dioleoyl-2-palmitoylglycerol

1,3-Di(cis-9-octadecenoyl)-2-hexadecanoylglycerol

1,3-Dioleo-2-palmitin

2-Palmito-1,3-diolein

CAS: 1716-07-0

Catalog #: RP2327, 25mg


Formula: $C_{55}H_{102}O_6$

MW : 859.39

Appearance: Clear Colourless Oil

Solubility: Chloroform (Slightly),
Methanol (Slightly)

Storage: Store at -20°C



References

- Nakaya, H. *et al.*: Biotechnology Techniques, 12, 881 (1998)
- Chen, M. *et al.*: J. Am. Oil Chemists' Society, 91, 525 (2004)

Ordering information

Catalog size quantities and prices may be found at <http://www.interchim.com>.

Please inquire for higher quantities (availability, shipment conditions).

Please contact InterBioTech – Interchim for any other information

Hotline : +33(0)4 70 03 73 06 – Interbiotech@interchim.com

Disclaimer : Materials are sold **for research use only**, and are not intended for food, drug, household, or cosmetic use. Interchim® is not liable for any damage resulting from handling or contact with this product

