

FT-U0786B



## RNase Removal Spray

*Removal of RNase & DNA contamination*

### Product Description

**Name :** **RNase Removal Spray**  
 RNase Removal Spray Solution for molecular biology.

**Catalog Number :** U0786B, 1L

**Storage:** Room temperature

### Introduction

This handy 1L spray solution removes RNase and DNA from laboratory surfaces and is effective against enzyme contamination. This will decontaminate bench tops, instruments, pipettors, labware, working well on glass or plastic surfaces.

Available as a ready to use solution simply spray on the area contaminated and wipe away from the surface using distilled water.

### Ordering information

Catalog size quantities and prices may be found at <http://www.interchim.com>. Please inquire for higher quantities (availability, shipment conditions).

For any information, please ask : Uptima / Interchim; Hotline : +33(0)4 70 03 73 06

[Order on-line](#) or [Contact](#) your local distributor

**Disclaimer :** Materials from Uptima are sold **for research use only**, and are not intended for food, drug, household, or cosmetic uses. Uptima is not liable for any damage resulting from handling or contact with this product.