

- Prevention is better than cure - (U)HPLC Guard Columns

Protection for 1.7 up to 5 μ m totally porous silicas

Maximum pressure : 900 bar # Low dead volume # High performance

Guard Holder
P/N : AGHP-5



Guard cartridge



Uptisphere® HP Guard cartridges	Reverse Phase	Hilic Mode	Normal Phase
1.7 μ m - 5 x 2.1mm - 3u	UP-RP-1.7-005/021	---	---
1.7 μ m - 5 x 4.0mm - 3u	UP-RP-1.7-005/046	---	---
2.2 μ m - 5 x 2.1mm - 3u	UP-RP-2.2-005/021	UP-HILIC-2.2-005/021	UP-NP-2.2-005/021
2.2 μ m - 5 x 4.0mm - 3u	UP-RP-2.2-005/046	UP-HILIC-2.2-005/046	UP-NP-2.2-005/046
3 μ m - 5 x 2.1mm - 3u	UP-RP-3-005/021	UP-HILIC-3-005/021	UP-NP-3-005/021
3 μ m - 5 x 4.0mm - 3u	UP-RP-3-005/046	UP-HILIC-3-005/046	UP-NP-3-005/046
5 μ m - 5 x 2.1mm - 3u	UP-RP-5-005/021	UP-HILIC-5-005/021	UP-NP-5-005/021
5 μ m - 5 x 4.0mm - 3u	UP-RP-5-005/046	UP-HILIC-5-005/046	UP-NP-5-005/046

ANALYTICAL SCIENCES ITM-20160119-C/P

 **interchim**®

Hotline +33 4 70 03 73 09 • interchrom@interchim.com


Suivez-nous sur :



- Prevention is better than cure - Core Shell Guard Columns

Protection for 2.6 - 2.7 μ m Core Shell Silica

Maximum pressure : 1000 bar # Low dead volume # High performance



Uptisphere [®] CS-HP Guard cartridges	Reverse Phase	Hilic Mode
5 x 2.1mm - 3 μ	UE-RP-2.6-005/021	UE-HILIC-2.6-005/021
5 x 3.0mm - 3 μ	UE-RP-2.6-005/030	UE-HILIC-2.6-005/030
5 x 4.6mm - 3 μ	UE-RP-2.6-005/046	UE-HILIC-2.6-005/046

ANALYTICAL SCIENCES ITM-20160119-C/P

 **interchim**[®]

Hotline +33 4 70 03 73 09 • interchrom@interchim.com

Suivez-nous sur :

