

FT-AXFCRA



## UptiLight Clean WB HRP Chemiluminescent Substrate

### Product Description

<b>Name :</b>	<b>UptiLight™ Clean WB HRP Chemiluminescent Substrate</b>
<b>Catalog Number :</b>	AXFCRA, 2 x 30 ml (for 600cm <sup>2</sup> membrane) AXFCRB, 2 x 60 ml (for 1200cm <sup>2</sup> membrane) AXFCRA, 2 x 125 ml (for 2500cm <sup>2</sup> membrane)
<b>Storage:</b>	+2-8 °C (12 months)

### Description


#### Applications

- Western Blotting
- Dot Blot Assay

#### Specifications

- Signal intensity: Standard
- Compatible with Nitrocellulose & PVDF membrane
- Detection of horseradish peroxidase (HRP) enzyme activity
- Color: colorless

Detection with UptiLight™ Clean WB HRP  
Chemiluminescent Substrate



### Ordering information

Catalog size quantities and prices may be found at <http://www.interchim.com>.  
Please inquire for higher quantities (availability, shipment conditions).

For any information, please ask : Uptima / Interchim; Hotline : +33(0)4 70 03 73 06

[Order on-line](#) or [Contact](#) your local distributor

**Disclaimer :** Materials from Uptima are sold **for research use only**, and are not intended for food, drug, household, or cosmetic uses.  
Uptima is not liable for any damage resulting from handling or contact with this product.