

Whatman®



MONTLUCON
211 bis av. J.F. Kennedy - B.P. N° 1140
03103 Montluçon Cedex - France
☎(33) (0)4 70 03 88 55 (+)
Fax (33) (0)4 70 03 82 60
interbiotech@interchim.com




3MM Chr Blotting Papers

Whatman 3MM Chr Paper is the most widely used blotting paper. This acceptance and usage reflect the high quality, purity, and consistency that are relied upon by researchers doing Southern, Northern, and Western transfers. 3MM Chr Paper is now available in the most commonly used sizes.

FEATURES

- **Pure Cellulose:** Whatman 3MM Chr Paper is made entirely from the highest quality cotton linters. There are no additives of any kind ever used. This ensures that no contamination will occur during the transfer steps, as can happen when using inferior papers.
- **Special Chromatographic Design:** Whatman 3MM Chr Papers are manufactured and tested specially for chromatographic techniques. This ensures the wicking capability and uniformity of capillary action that is so important in obtaining clean and even transfers during blotting. Just because a paper is thick doesn't guarantee it is a good blotting material.
- **Consistent Quality:** Whatman 3MM Chr Paper was used in the earliest blotting procedures. Manufacturing to specifications for chromatography paper have made 3MM Chr the consistent standard in the industry.
- **Convenient Sizes:** Whatman 3MM Chr Papers are now available in sheets precisely cut to the most popular gel and transfer membrane sizes. This allows "out-of-the-box" usage and eliminates sheet-to-sheet variations.
- **Easy Ordering Procedure:** Whatman 3MM Chr Papers are available through our extensive network of national and regional distributors. You can purchase any of these products in the same manner as your other lab supplies and receive the superior service to which you are accustomed.

3MM Chr Paper Typical Blotting Apparatus



Whatman Quality

Whatman is a global leader in separations technology and is known in the scientific community for providing Innovative products and solutions. Our instinct for simplification accelerates the rate of discovery, reduces costs and saves time. For more information, visit www.whatman.com

Whatman® and Protran® are trademarks of the Whatman Group.

North America

Whatman Inc.
200 Park Avenue, Suite 210
Florham Park, NJ 07932 USA
WHATMAN (US and Canada)

Europe



MONTLUCON
211 bis av. J.F. Kennedy - B.P. N° 1140
03103 Montluçon Cedex - France
☎(33) (0)4 70 03 88 55 (+)
Fax (33) (0)4 70 03 82 60
interbiotech@interchim.com

3MM Chr Pure Cellulose Grade Sheets

Catalog Number	Size	Quality/Packs
3030-6185	11 x 14 cm	100
3030-6366	11 x 20 cm	100
3030-6132	12 x 14 cm	100
3030-6463	12 x 21 cm	100
3030-153	15 x 17.5 cm	100
3030-6188	15 x 20 cm	100
3030-221	18 x 34 cm	100
3030-861	20 x 20 cm	100
3030-6461	26 x 41 cm	100
3030-347	35 x 43 cm	100
3030-392	35 x 45 cm	100
3030-917	46 x 57 cm	100
3030-931	58 x 68 cm	100
3030-6186	2.5" x 9"	100
3030-6481	3.25" x 5"	100
3030-6189	4" x 5.25"	100
3030-6187	6" x 8"	100
3030-866	8" x 10"	100

*To place an order for Whatman Protran® Nitrocellulose Blotting Membrane, please call your local Whatman distributor.