

Cat.-No.	Size
AB-104S	1 ml
AB-104L	5 ml

Description:

Contains 100 mg/ml LEXSY Hygro antibiotic.
Sterile. Based on antibiotic Hygromycin B TM.

Storage conditions:

Store at -20°C.
Stable for 18 month.

Selection:

For selection of recombinant *Leishmania* strains
LEXSY Hygro is added to growth medium to a
final concentration of 100 µg/ml.

Find more information in the User's Guide for
Leishmania tarentolae Gene Expression Starter Kit
(download pdf from the LEXSY section of our
website).

FOR RESEARCH USE ONLY.

**NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC
OR THERAPEUTIC USES.**