



TRIS- EDTA 100X

Product Description

Name :	TRIS-EDTA Buffer 100X, Biotechnology Grade (TE BUFFER 100X)
Catalog Number :	IOY832, 500 ml
Description:	A single strength (1X) working solution contains 10 mM Tris, 1mM EDTA
Storage:	+4°C (L)

Tris-EDTA Buffer is commonly used as storage buffer for DNA and RNA samples. The Buffer protects DNA and RNA from degradation.

Ordering information

Catalog size quantities and prices may be found at <http://www.interchim.com>.

For any information, please ask : Uptima / Interchim; Hotline : +33(0)4 70 03 73 06

[Order on-line](#) or [Contact](#) your local distributor

Disclaimer : Materials from Uptima are sold **for research use only**, and are not intended for food, drug, household, or cosmetic uses. Uptima is not liable for any damage resulting from handling or contact with this product.

Rev.P01vb

Uptima™, powered by

P.1