

PENTAPHARM

# Pefachrome<sup>®</sup> ELA

**Application:** Highly sensitive chromogenic peptide substrate for leucocyte elastase. Determination of elastase activity for research, in-process and quality control.

**Formula:** MeOSuc-Ala-Ala-Pro-Val-pNA

**Principle:** MeOSuc-Ala-Ala-Pro-Val-pNA + E ==> MeOSuc-Ala-Ala-Pro-Val-OH + pNA + E  
E = elastase

**MW:** 590.6

**Storage:** May be used by the expiry date given on the label when stored unopened, protected from moisture, in the dark, 2-8°C. Avoid contamination of the reagents by micro-organisms.  
Shipment of product does not require cooling during the time of transportation.

**Package size:** Vial containing 25 mg  
Bulk [g]

**Code:** 090-08

FOR RESEARCH USE ONLY. NOT FOR HUMAN USE OR DRUG USE.

