


FT-JO3540



Sequencing Gel 19:1

Sequencing Gel 19:1 ratio

Product Description

Name :	6% Sequencing Gel 19:1 ratio
Catalog Number :	6% Manual Sequencing Gel Solution Bis-Acrylamide Ratio 19:1 (Ready to use formulation) JO3540, 500ml JO3541, 1L
Name :	8% Sequencing Gel 19:1 ratio
Catalog Number :	8% Manual Sequencing Gel Solution Bis-Acrylamide Ratio 19:1 (Ready to use formulation) 1Q6690, 1L 1Q6691, 500 ml
Name :	4.5% Sequencing Gel 19:1 ratio
Catalog Number :	4.5% Manual Sequencing Gel Solution Bis-Acrylamide Ratio 19:1 (Ready to use formulation) 1R0360, 500 ml 1R0361, 1 L

Storage: +2-8°C Minimum shelf life 6 months

Used in routine DNA sequencing work, ready to use formulation contains 7M Urea (Denaturant) and 1X TBE buffer. Each lot is DNase/ RNase tested. All stock solutions are made from ultra-pure components and both 0.45 and 0.2micron filtered.

Ordering information

Catalog size quantities and prices may be found at <http://www.interchim.com>. Please inquire for higher quantities (availability, shipment conditions).

For any information, please ask : Uptima / Interchim; Hotline : +33(0)4 70 03 73 06
[Order on-line](#) or [Contact](#) your local distributor

Disclaimer : Materials from Uptima are sold **for research use only**, and are not intended for food, drug, household, or cosmetic uses. Uptima is not liable for any damage resulting from handling or contact with this product.