

FT-JW7390




## DMSO anhydrous

*Organic solvent used for dissolving moisture-sensitive compounds*

### Product Description

<b>Name :</b>	<b>Dimethyl sulfoxide, anhydrous (DMSO)</b>
<b>Catalog Number :</b>	FP-JW7390, 10ml
<b>Molecular Formula:</b>	(CH <sub>3</sub> ) <sub>2</sub> SO
<b>MW :</b>	78.13
<b>Color &amp; Form:</b>	Colorless liquid
<b>Purity:</b>	≥ 99.9%
<b>Titration :</b>	0.005% H <sub>2</sub> O (maximum)
<b>Residual on Evaporation :</b>	0.0005% (maximum)
<b>Storage :</b>	Store desiccated at room temperature



### Ordering information

Catalog size quantities and prices may be found at <http://www.interchim.com>.  
Please inquire for higher quantities (availability, shipment conditions).

For any information, please ask : FluoProbes® / Interchim; Hotline : +33(0)4 70 03 73 06

**Disclaimer :** Materials from FluoProbes® are sold **for research use only**, and are not intended for food, drug, household, or cosmetic use. FluoProbes® is not liable for any damage resulting from handling or contact with this product.